



Behind Every Smile is a Story

- QUALITY - ЯКІСТЬ
- PRECISION - ТОЧНІСТЬ
- ACCURACY - ПРАВИЛЬНІСТЬ
- DURABILITY - ЗНОСОСТІЙКІСТЬ



## ORTHODONTIC INSTRUMENTS

- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente
- Инструменти ортодонтични
- Instruments Orthodontiques
- أدوات تقويم الأسنان



**AG5-1800**  
Height Gauge 0.18  
Aluminum



**AG5-1801**  
Height Gauge 0.22  
Aluminum



**AG5-1802**  
Height Gauge 0.18  
Stainless Steel



**AG5-1803**  
Height Gauge 0.22  
Stainless Steel



**AG5-1804**  
Moveable  
Height Gauge



**AG5-1805**  
Endo Finger  
Ruler Span  
Measure Scale



**AG5-1806**  
Endo Finger  
Ruler Span  
Measure Scale



**AG5-1807**  
Endo Finger  
Ruler Span  
Measure Scale



**AG5-1808**  
Positive  
Postioner



**AG5-1809**  
Strip Saw



**AG5-1810**  
Bone Bracket Gauge



**AG5-1811**  
Endo Finger Ruler  
Span Measure Scale



**AG5-1812**  
Endo Finger Ruler  
Span Measure Scale

## ORTHODONTIC INSTRUMENTS

- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente
- Инструменти ортодонтични
- Instruments Orthodontiques
- ادوات تقويم الأسنان



**AG5-1813**  
Measuring Loop  
angled, length  
175 mm,  
measuring range  
0 - 40 mm



**AG5-1814**  
Measuring Loop  
straight length  
175 mm,  
measuring range  
0 - 40 mm



**AG5-1815**  
Measuring loop, Castroviejo  
angled, length 175 mm,  
measuring range 0 - 40 mm,  
with window



**AG5-1816**  
Measuring Loop  
straight length 90 mm,  
measuring range  
0 - 40 mm,  
**AG5-1817**  
Castroviejo Calipers  
20mm  
90mm



**AG5-1818**  
Measuring Loop  
angled, length 90 mm,  
measuring range  
0 - 40 mm,  
**AG5-1819**  
Castroviejo Calipers  
20mm, 45°  
90mm



**AG5-1820**  
Castroviejo Calipers  
20mm  
90mm



**AG5-1821**  
Fig. 1  
Iwanson  
spring caliper  
For metal and  
porcelain.



**AG5-1822**  
Fig. 2  
Iwanson  
spring caliper  
For wax.

## ORTHODONTIC INSTRUMENTS


- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente
- Инструменти ортодонтични
- Instruments Orthodontiques
- ابوات تقويم الأسنان



**AG5-1823**  
Boley Gauge



**AG5-1824**  
Vernier Gauge



**AG5-1825**  
Aluminium Caliper  
up to 200mm



**AG5-1826**  
Boley Gauge



**AG5-1827**  
Delux Boley  
Gauge



**AG5-1828**  
Beerendonk  
130mm  
Mini Gauge



**AG5-1829**  
Compasso  
De Willis



**AG5-1830**  
Dentaurum  
Ruler



**AG5-1831**  
Ruler Wedge



**AG5-1832**  
Ruler

## ORTHODONTIC INSTRUMENTS

- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente
- Инструменти ортодонтични
- Instruments Orthodontiques
- ادوات تقويم الأسنان

### Bracket Holder Tweezers



**AG5-1833**  
Deluxe Bracket  
Holder Tweezers



**AG5-1834**  
Deluxe Bracket  
Holder Tweezers



**AG5-1835**  
Direct Bond Bracket  
Holder Tweezers



**AG5-1836**  
Bracket Tweezers  
Slot Aligner



**AG5-1837**  
Direct Bond  
Bracket Holder  
Tweezers, Small Tip



**AG5-1838**  
Exact Bracket  
Holder Tweezers



**AG5-1839**  
Bracket Placing  
Tweezers



**AG5-1840**  
Easy Access  
Tweezers



**AG5-1841**  
Bracket Placer  
Tweezers  
(Extra Fine Tip)



**AG5-1842**  
Deluxe Bracket  
Holder Tweezers  
135mm



**AG5-1843**  
Easy Eccess  
Tweezers



**AG5-1844**  
Slipper  
Tweezers



**AG5-1845**  
Direct Bonding  
Bracket Holders  
135mm



**AG5-1846**  
Direct Bonding  
Bracket Holders  
136mm



**AG5-1847**  
Direct Bonding  
Bracket Holders  
135mm



**AG5-1848**  
Deluxe Bracket  
Holder



**AG5-1849**  
Direct Bonding  
Bracket Holders  
140mm

## ORTHODONTIC INSTRUMENTS

- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente
- Инструменти ортодонтични
- Instruments Orthodontiques
- ابوات تقويم الأسنان



**AG5-1850**  
Band seater  
auto cleavable



**AG5-1851**  
Band seater  
auto cleavable  
Triangle Shape



**AG5-1852**  
Band seater  
auto cleavable  
Square Shape



**AG5-1853**  
Band seater  
auto cleavable  
Round Shape



**AG5-1854**  
Mershon Band  
Pusher  
140mm



**AG5-1855**  
Band seating &  
removing  
instruments



**AG5-1856**  
Ortho Light double  
ended  
ligature director  
scaler



**AG5-1857**  
Ortho Light Band  
Pusher Scaler



**AG5-1858**  
Band Seating &  
Removing  
Silicon Grip



**AG5-1859**  
Ligature Instrument  
Silicon Grip



**AG5-1860**  
Elastic Instrument  
Silicon Grip



**AG5-1861**  
Distal Bender  
Silicon Grip



**AG5-1862**  
Distal Bender  
Silicon Grip



**AG5-1863**  
Distal Bender  
Silicon Grip

## ORTHODONTIC INSTRUMENTS

- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente
- Инструменти ортодонтични
- Instruments Orthodontiques
- ادوات تقويم الأسنان



**AG5-1864**  
Ortho light scaler double ended



**AG5-1865**  
Ortho light scaler double ended



**AG5-1866**  
Ortho light scaler double ended



**AG5-1867**  
Twiel-on elastic module instrument



**AG5-1868**  
Twiel-on elastic module instrument



**AG5-1869**  
Schure Band Seater



**AG5-1870**  
Ortho light scaler



**AG5-1871**  
Ortho light scaler



**AG5-1872**  
Ortho light scaler heavy



**AG5-1873**  
Ortho instrument



**AG5-1874**  
Ortho instrument

## ORTHODONTIC INSTRUMENTS

- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente
- Инструменти ортодонтични
- Instruments Orthodontiques
- ابوات تقويم الأسنان



**AG5-1875**  
Ortho Light double ended  
ligature director



**AG5-1876**  
Ortho Light double ended ligature scaler director tivate coil spring  
are remove modules with one instruments



**AG5-1877**  
Ortho band pusher scaler



**AG5-1878**  
Ortho ligature director single ended



**AG5-1879**  
Ortho Light double ended ligature  
director scaler



**AG5-1880**  
Ortho Light double ended ligature scaler/director  
tivate coil spring are remove modules with one instruments



**AG5-1881**  
Ortho Light double ended ligature director



**AG5-1882**  
Ortholite Ligature Director Scaler



**AG5-1883**  
Ortho light ligature director double ended



**AG5-1884**  
Double ended cement spatula



**AG5-1885**  
Ortho light scaler single ended



## PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمشة لتقويم الأسنان برومبئیس

### BOX JOINT PLIERS

**AG5-1886**  
Universal Dolphin Style  
Wire Bending & cutting pliers



**AG5-1892**  
Bracket Removing Plier curved  
Bracket Removing Plier curved for  
easy removal of brackets in the  
posterior region. Length 130mm



**AG5-1887 Band Removing Plier**  
The curved blade of the band  
removing plier can strip easily  
and gently remove directly to  
the tooth. Length 130mm



**AG5-1893**  
Young Wire Bending &  
Loop Forming Pliers



**AG5-1888**  
Tweed Loop  
forming Pliers



**AG5-1894**  
Wire Bending &  
Loop Forming Pliers



**AG5-1889**  
How Plier curved  
Universally usable tool for  
the multiband technology for  
wires up to .028".  
Length 130mm



**AG5-1895**  
Adams, small wire bending  
and loop forming Pliers



**AG5-1890 How Plier**  
Universally usable tool for the  
multiband technology for wires  
up to .028".  
Length 135mm



**AG5-1896**  
Bracket Removing Plier Round Box  
Bracket removing plier for comfortable  
and quick removal of the brackets.  
Length 135mm



**AG5-1891**  
Nance Loop  
forming Pliers



**AG5-1897**  
Weingart Plier medium  
Weingart plier that are ideal due to  
their curved plier tip for setting sheets and  
is capable of bending the ends.  
Length 145mm




## PLIERS FOR ORTHODONTICS & PROSTHETICS


- Щипцы для ортодонтии и протезирования
- Щипці для ортодонції та протезування
- Pincers pour l'orthodontie et la prothèse
- Alicates para ortodondia y prótesis
- Zangen für Kieferorthopädie und Prothetik
- كمشة لتقويم الأسنان بروتستیس

### BOX JOINT PLIERS

**AG5-1898 Adams Plier Round Box**  
Adams pliers for bending wire arches and bends right angles for wires up to .022", .025 tip to beak. Length 125mm



**AG5-1904 Adams Plier**  
Adams pliers for bending wire arches and bends right angles for wires up to .022", .025 tip to beak. Length 125mm



**AG5-1899 Falcon Beak Contouring Pliers.**



**AG5-1905 Tweed Loop Formin Pliers**



**AG5-1900**  
wires up to 0.020" / 0.7mm.



**AG5-1906 Bird Beak Plier Standard**  
Bird Beak Wire bending plier without notches, which is suitable as a universal utility tool. Length 125mm



**AG5-1901 Johnson Contouring Pliers**



**AG5-1907 Mini Ligature cutters**



**AG5-1902 Aderer Large wire bending and loop forming pliers**



**AG5-1908 Small Fluch-cut distal and cutting Pliers**  
with safety hold fitted with un-breakable, corrosion resistant and durable stainless steel tips



**AG5-1903 Tweed / angle rectangular arch forming & Torquing pliers**



**AG5-1909 Large Fluch-cut distal and cutting Pliers**  
with safety hold fitted with un-breakable, corrosion resistant and durable stainless steel tips



## PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمشة لتقويم الأسنان برومینیس

### BOX JOINT PLIERS

**AG5-1910**  
**Tweed / Angel**  
1.78mm (.070") wide tips arch  
forming and torquing pliers



**AG5-1916**  
Contouring Plier



**AG5-1911**  
**Tweed / Angel**  
1.78mm (.070") wide tips arch  
forming and torquing pliers



**AG5-1917**  
Bands & Crowns Plier  
in order to skirt



**AG5-1912**  
Lingual / Palatal arch  
forming pliers



**AG5-1918**  
Ligature Forming Pliers



**AG5-1913**  
Aderer  
Facebow loop forming Pliers



**AG5-1919**  
for double back and triple back  
bands in 0.30" / 0.76mm & 0.36"  
/ 0.91mm lingual arch wires



**AG5-1914**  
Waldsavhs wire bending  
and cutting pliers



**AG5-1920**  
Bird Beak Wire Bending  
& cutting pliers



**AG5-1915 Aderer Plier**  
For the manufacturing of appliances and  
activation of wires for wires up to .030" .fine  
tip and up to .036" mean tip.  
Length 130mm



**AG5-1921**  
Weingart Plier medium  
Weingart plier that are ideal due to  
their curved tip for setting sheets and  
is capable of bending the ends.  
Length 145mm



## PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincas pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كاشطة لتقويم الأسنان بروسیتیس

### L-KEY JOINT PLIERS



**AG5-1922**  
Distal end cutter with safety hold, mini head, long handles.  
fitted with un-breakable, corrosion resistant and durable stainless steel tips

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.46 mm x 0.64 mm (.018" x .025")
- ⊗ 0.38 mm to 0.53 mm (.0155" to .0215")
- ⊗ Max. 0.53 mm x 0.64 mm (.021" x .025")



**AG5-1923**  
Flush-cut distal end cutter with safety hold fitted with un-breakable, corrosion resistant and durable stainless steel tips

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.46 mm x 0.64 mm (.018" to .025")



**AG5-1924**  
Soft wire and ligature cutter, Mini fitted with un-breakable, corrosion resistant and durable stainless steel tips

- ⌋ .25 mm to .38 mm (.010" to .015")



**AG5-1925**  
Soft wire and ligature cutter, standard fitted with un-breakable, corrosion resistant and durable stainless steel tips

- ⌋ .25 mm to .38 mm (.010" to .015")



**AG5-1926**  
Soft wire and ligature cutter, 15° angled fitted with un-breakable, corrosion resistant and durable stainless steel tips

- ⌋ .25 mm to .38 mm (.010" to .015")



**AG5-1927**  
Soft wire and ligature cutter, 45° angled fitted with un-breakable, corrosion resistant and durable stainless steel tips

- ⌋ .25 mm to .38 mm (.010" to .015")



**AG5-1928**  
Hard wire cutter, straight Fitted with un-breakable, corrosion resistant and durable stainless steel tips

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



**AG5-1929**  
Hard wire cutter, 15° angled Fitted with un-breakable, corrosion resistant and durable stainless steel tips

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")
















**AG5-1930**  
Hard wire cutter, 15° angled Non-autoclavable

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")

## PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمشة لتقويم الأسنان بروستیس

### L-KEY JOINT PLIERS

	 <p><b>AG5-1931</b> Band removing pliers</p>		 <p><b>AG5-1936</b> Weingart-mini Utility Pliers</p> <ul style="list-style-type: none"> <li>● 0.30 mm to 0.50 mm (.012" to .020")</li> <li>■ Max 0.56 mm x 0.71 mm (.022" x .028")</li> </ul>
	 <p><b>AG5-1932</b> Bracket removing pliers angled</p>		<p><b>AG5-1937</b> Weingart Utility Pliers</p> <ul style="list-style-type: none"> <li>● 0.30 mm to 0.50 mm (.012" to .020")</li> <li>■ Max 0.56 mm x 0.71 mm (.022" x .028")</li> </ul>
	 <p><b>AG5-1933</b> Bracket removing pliers straight</p>		 <p><b>AG5-1938</b> How Straight Utility Pliers</p> <ul style="list-style-type: none"> <li>● 0.30 mm to 0.50 mm (.012" to .020")</li> <li>■ Max 0.56 mm x 0.71 mm (.022" x .028")</li> </ul>
	 <p><b>AG5-1934</b> Bracket removing pliers with plastic tip</p>		<p><b>AG5-1939</b> Angled Utility Pliers</p> <ul style="list-style-type: none"> <li>● 0.30 mm to 0.50 mm (.012" to .020")</li> <li>■ Max 0.56 mm x 0.71 mm (.022" x .028")</li> </ul>
	 <p><b>AG5-1935</b> Adhesive removing pliers</p>		 <p><b>AG5-1940</b> Wire bending and loop forming pliers Angled (Bird beak)</p> <ul style="list-style-type: none"> <li>● 0.30 mm to 0.50 mm (.012" to .020")</li> <li>■ Max 0.56 mm x 0.71 mm (.022" x .028")</li> </ul>

## PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Пинцес pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كاشطة لتقويم الأسنان بروتستیس

### L-KEY JOINT PLIERS



**AG5-1941**  
Multi-purpose pliers, bird beak,  
De la rosa and cutter in one

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.46 mm x 0.56 mm (.018" x .022")



**AG5-1946**  
Stop pliers V-Bebd - 1mm

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.43 mm x 0.56 mm (.017" x .028")



**AG5-1942**  
Cinch Back / tie back  
Distal end bending pliers

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



**AG5-1947**  
Lingual / palatal arch forming plies

- 0.76 mm (.030") / 0.91 mm (.036")



**AG5-1943**  
Hammer head  
Distal end bending pliers

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



**AG5-1948**  
Torquing pliers set



**AG5-1949**  
Tweed / angel  
Rectangular arch forming and  
torquing pliers

- 0.30 mm to 0.71 mm (.012" to .028")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



**AG5-1944**  
Nance  
Loop forming pliers

- 0.36 mm to 0.71 mm (.014" to .028")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



**AG5-1950**  
Surgical ball hook crimping pliers,  
angled



**AG5-1945**  
Detailing step pliers / Bayonet  
bending pliers - 1mm

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



**AG5-1951**  
Surgical ball hook crimping  
pliers, straight



## PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمشة لتقويم الأسنان برومینیس


### L-KEY JOINT PLIERS





**AG5-1952**  
Aderer Clamp three jaw wire bending plier for contouring and bending wire up to .30"/0.76mm.



**AG5-1958**  
Falcon beak contouring plier.




**AG5-1953**  
Bird beak plier for bending round, square and rectangle wire up to .30"/0.76mm.



**AG5-1959**  
For placement and removal of arch wires, pins and other auxiliaries.




**AG5-1954**  
Bird beak plier for light wire bending round, square and rectangular up to .016"/0.41mm.



**AG5-1960**  
For placement and removal of arch wires, pins and other auxiliaries in angle design



**AG5-1955**  
Double Rounded jaw bending pliers, hooks and wire up to .020"/0.51mm.




**AG5-1961**  
contouring Plier.



**AG5-1962**  
Bending plier. In order to skirt molares bands. With grooved. Recommended wire. 016", 018" and 022" 0.41mm, 0.46mm and 0.56mm



**AG5-1956**  
Adams pliers for hooks and bending wire up to .030"/0.76mm.



**AG5-1963**  
Loop forming in four length 3,4,5 & 6 mm wire capacity 0.22", 56mm.



**AG5-1957**  
Johnson pliers for crown band contouring



**AG5-1964**  
Arch bending plier design for convenient handling of square or rectangular wire, .020" / 0.51mm.



# PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincas pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodondia y prótesis
- كمشة لتقويم الأسنان بروستیتیس


## LAP & L-KEY JOINT PLIERS



**AG5-1965**  
In order to skirt bands and crowns.



**AG5-1966**  
Wires up to .020" / 0.7mm.



**AG5-1967**  
Round and concave plier for loop forming and banding wire up to .20" / .051mm.




**AG5-1968**  
Bending for round square or rectangular wire up to 0.30"/0.76mm



**AG5-1969**  
for double back and triple back bands in 0.30" / .76mm & 0.36" / .91mm lingual arch wires.



**AG5-1970**  
Detailing step plier double-sided offset beak permits bayonet bends. Creates right hand or left hand step-up or stepdown at 1/4, 3/4, or 1mm increments.



**AG5-1971**  
"V" bend stop special design makes 1mm "v" bends to shorten archwires or provide a positive stop. Excellent for placing stops in titanium wires. For round or rectangular wires up to .025"



**AG5-1972**  
Angle how plier efficient utility plier with serrated beaks for a firm grip. 45 angled tips to access hard to reach areas.



**AG5-1973**  
Weingardt plier serrated tips hold wires firmly at a convenient working angle. Tapered beaks fit easily between brackets and hard to reach areas. Our most popular utility plier.



**AG5-1974**  
Pin & Ligature Cutters TC.



**AG5-1975**  
Mini Pin & Ligature Cutter TC.



**AG5-1976**  
Standard Ligature Cutter TC.



**AG5-1977**  
V Bend Stop Plier TC.



**AG5-1978**  
Angled Direct Bond Bracket Remover TC.



**AG5-1979**  
Cutter Str.



Distal Cutter 14cm  
**AG5-1980** 2.5 mm  
**AG5-1981** 3.0 mm  
**AG5-1982** 3.0 mm



**AG5-1983**  
Arrow Clasp Forming Pliers 12.0cm



**AG5-1984**  
Schwarz 12.5cm



## DISTAL END, WIRE CUTTERS & BRACKET REMOVER PLIERS

**AG5-1985**  
Distal end Cutter, Safety Hold Plier  
(Ideal for Cutting Wire) TC.



**AG5-1986**  
Hook Crimping Plier Angled  
Straight TC.



**AG5-1987**  
3 Prong Clasp Adjustable  
Plier TC.



**AG5-1988**  
Utility Plier Weingart TC.



**AG5-1989**  
Direct Bracket Remover  
TC.



**AG5-1990**  
Straight Direct Bond Bracket  
Remover TC.



**AG5-1991**  
Pin Holding Plier TC.



**AG5-1992**  
Hard wire cutters



**AG5-1993**  
Elastic Separating Plier  
Round Handle



**AG5-1994**  
Elastic  
Separating Plier



**AG5-1995 Str.**  
**AG5-1996 Cvd.**  
Smaha  
Elastomeric  
ligature  
placing forceps  
140 mm



**AG5-1997**  
Elastomeric  
ligature placing  
forceps with  
groove  
Falcon 140 mm



**AG5-1998**  
Elastomeric  
ligature placing  
forceps with  
groove  
Falcon 140 mm



## PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincas pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كماشة لتقويم الأسنان بروسیتیس

### CUTTER, UTILITY PLIERS & WIRE FORMING PLIERS




**AG5-1999**  
Weingart Pliers



**AG5-2004**  
Adhesive  
Removing Pliers



**AG5-2000**  
Slim Weingart  
Pliers



**AG5-2005**  
Posterior Band Removing  
Pliers - Long



**AG5-2001**  
How Pliers  
Straight



**AG5-2006**  
Posterior Band Removing  
Pliers - Short



**AG5-2002**  
How Pliers Offset




**AG5-2007**  
Bird Beak Pliers



**AG5-2003**  
Band seating Pliers



Replacement tip and  
hex key included.



**AG5-2008**  
Omega Loop  
Forming Pliers

## PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمشة لتقويم الأسنان بروستیس

### CUTTERS



**AG5-2009**  
Hard  
Wire Cutter 15



**AG5-2013**  
Mini Pin &  
Ligature Cutter



**AG5-2010**  
Distal End Cutter  
Flush Cut



**AG5-2014**  
Hard Wire Cutter  
Straight



**AG5-2011**  
Universal Cut & Hold  
Distal End Cutter



**AG5-2015**  
Micro-Mini Pin &  
Ligature Cutter



**AG5-2012**  
Pin & Ligature  
Cutter



**AG5-2016**  
Cutter

## PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Щипці для ортодонції та протезування
- Pincas pour l'orthodontie et la prothèse
- Alicates para ortodoncia y prótesis
- Zangen für Kieferorthopädie und Prothetik
- كماشة لتقويم الأسنان بروتستیس



**AG5-2017**  
Three Jaw Pliers



**AG5-2021**  
Jarabak Pliers



**AG5-2018**  
Light Wore Bird Beak  
Pliers with Groove



**AG5-2022**  
Hollow Chop  
Contouring Pliers



**AG5-2019**  
Tweed loop  
Forming Pliers



**AG5-2023**  
Lingual Arch  
Forming Pliers



**AG5-2020**  
Arch Bending Pliers



**AG5-2024**  
Rectangular Arch  
Bending Pliers

## ORTHODONTICS PLIERS STANDS

- Ортодонтия Щипцы Подставки
- Supports de pinces d'orthodontie
- Kieferorthopädie Zangenständer
- Ортодонтия Щипці Підставки
- KFO Zangen Ständer
- تقويم الأسنان كمشاة تقف



**AG5-2025**  
Plier Stand  
Stainless Steel.



**AG5-2026**  
Plier Stand  
Stainless Steel.



**AG5-2027**  
Pliers Stand  
Stainless Rod  
(Aluminum Anodized)  
Stainless Steel



**AG5-2028**  
Plier Stand  
Stainless Steel.



**AG5-2029**  
Plier Stand  
Stainless Steel.



**AG5-2030**  
Plier Stand  
Plastic

ALLIED GROUP  
of Industries