


Our Scalers

Compare our scalers and curettes with the best in the industry and you will immediately recognize we are equal in every way but one... price. Note the deep precision knurling for better tactile feel. The same high-grade stainless steel tips, the same light weight (under 18 grams) for less hand fatigue, the same professional craftsmanship, and surpassing the highest industry standards. of course, any manufacturer can make these claims, but we guarantee and certify ours.



- High-polish mirror finish tips.
- Manufactured within the guidelines of the FDA quality assurance system for good manufacturing practices.
- Highest grade German stainless steel, razor sharp, durable tips. Rockwell hardened to 60 RHc.
- 27 of the most popular designs from gracesys to McCalls to columbias, and many european patterns.
- All instruments are guaranteed against defects in material and workmanship.

- All instruments bear the ce mark, emblematic of the european community's highest standards.
- Customer satisfaction guarantee, too!
- Less than 20 grams ... Lightweight for less hand fatigue.
- Manufactured to the most exacting standards so the 50th instrument you buy is exactly like your first.
- You decide ... Why pay more for the same high-quality instruments?

DIAGNOSTIC, PERIODONTAL & RESTORATIVE INSTRUMENTS

VALUED BUYERS

Compare our scalers and curettes with the best in the industry and you will immediately recognize we are equal in every way but one... price. Note the deep precision knurling for better tactile feel. The same high-grade stainless steel tips, the same light weight (under 20 grams) for less hand fatigue, the same professional craftsmanship, and surpassing the highest industry standards. Of course, any manufacturer can make these claims, but we guarantee and certify ours. We regularly manufacturing below articles in showing handles & pattern A to L for better grip & perfection. We are also fully equipped & capable to develop these instruments as per your own specific requirements & demands in materials, sizes, shapes and patterns.



A-Hollow Handle Ø 10 mm



B-Hollow Handle Ø 10 mm



C-Hollow Handle Ø 10 mm



D-Hollow Handle Ø 10 mm



E-Hollow Handle Ø 10 mm



F-Hollow Handle Ø 10 mm



G-Hollow Handle Ø 10 mm



H-Hollow Handle Ø 6.4 mm



I-Hollow Handle Ø 8 mm



J-Hollow Handle Ø 10 mm



K-Soild Handle Ø 6 mm



L-Soild Handle Ø 6 mm

MOUTH MIRRORS

- зеркала
 - Espejos de la boca
 - Mundspiegel
- дзеркала
 - Miroirs à la bouche
 - مرايا الفم

Plain



AG3-600 Fig. 3 Ø20mm
AG3-601 Fig. 4 Ø22mm
AG3-602 Fig. 5 Ø24mm
AG3-603 Fig. 6 Ø26mm

Cone-Socket Plain



AG3-604 Fig. 3 Ø20mm
AG3-605 Fig. 4 Ø22mm
AG3-606 Fig. 5 Ø24mm

Magnifying



AG3-607 Fig. 3 Ø20 mm
AG3-608 Fig. 4 Ø22 mm
AG3-609 Fig. 5 Ø24 mm

Cone-Socket Magnifying



AG3-610 Fig. 3 Ø20mm
AG3-611 Fig. 4 Ø22mm
AG3-612 Fig. 5 Ø24mm

Rhodium



AG3-613 Fig. 3 Ø20mm
AG3-614 Fig. 4 Ø22mm
AG3-615 Fig. 5 Ø24 mm

Cone-Socket Rhodium



AG3-616 Fig. 3 Ø20mm
AG3-617 Fig. 4 Ø22mm
AG3-618 Fig. 5 Ø24mm

Front Surface



AG3-619 Fig. 3 Ø20mm
AG3-620 Fig. 4 Ø22mm
AG3-621 Fig. 5 Ø24mm
AG3-622 Fig. 6 Ø26mm

Cone-Socket Front Surface



AG3-623 Fig. 3 Ø20mm
AG3-624 Fig. 4 Ø22mm
AG3-625 Fig. 5 Ø24mm

HANDLES FOR MOUTH MIRRORS

- Ручки для зеркал
- Manches pour miroirs à bouche
- Mundspiegelgriffe
- Ручки для дзеркал
- Mangos para espejos de boca
- مقابض لمرآة الفم



AG3-626
Hollow Mirror
Handle Round
Ø 10 mm



AG3-627
Hollow Mirror
Hexagonal
Handle
Ø 8 mm



AG3-628
Hollow
Ø 8 mm



AG3-629
Hollow Mirror
Hexagonal
Handle
Ø 10 mm



AG3-630
Hollow Mirror
Handle
Ø 8 mm



AG3-631
Hollow Mirror
Handle
Ø 10 mm



AG3-632
Mirror Handle
Solid
Ø 6 mm



AG3-633
Mirror Handle
Solid
Ø 6 mm



AG3-634
Mirror Handle
Solid Octagonal
Ø 6 mm



AG3-635
Mirror Handle
Stem
with Rule
35mm



AG3-636
Mirror Handle
Solid
Stem with
Rule 35mm



AG3-637
Cone-Socket
Solid Octagonal
Ø 6 mm



AG3-638
Cone-socket
Solid
Ø 6 mm



AG3-639
Cone-socket
Hexagonal
Handle
Ø 6 mm

Cone-Socket Mirror Handle

EXPLORERS

- Зонды
 - Sondes
 - Zahnsonden
- Зонди
 - Sondas
 - المستكشرون

Examination kit Set of 3 Pieces



AG3-640
Probe

AG3-641
Mirror

AG3-642
College Forcep
15cm



Examination Set of 4 Pieces



AG3-643
Exp,23/17-A


AG3-644
Excavators

AG3-645
Tweezers

AG3-646
Mouth Mirror

AG3-647

First Aid Examination Set of 3 Pieces



AG3-648

First Aid Examination Set of 3 Pieces



AG3-649

First Aid Examination Set of 3 Pieces



AG3-650

First Aid Examination Set of 3 Pieces



TWEEZERS

- Пинцеты
- Précettes - Brucelles
- Zahnpinzetten

- Пинцети
- Pinzas para algodón
- ملاقط



AG3-651
Surgical
11.5cm



AG3-652
London College
15cm



AG3-653
Semken Taylor
Str.



AG3-654
Semken Taylor
Cvd.



AG3-655
Meriam
16cm



AG3-656
Dressing
14cm



AG3-657
Semken Taylor
1x2 Straight



AG3-658
Semken Taylor
1x2 Curved



AG3-659
Gilles
15cm 1x2



AG3-660
14cm 1x2



AG3-661
Perry
13cm



AG3-662
Atraumatic
15cm Cvd



AG3-663
Atraumatic
15cm Str.



AG3-664
College with Lock
15cm Cvd.



AG3-665
Gerald
17.5cm



AG3-666
with lock
15cm



AG3-667
Atraumatic
17.5cm

TWEEZERS

- Пинцеты
- Précelles - Brucelles
- Zahnpinzetten

- Пинцети
- Pinzas para algodón
- ملاقط



AG3-668
Adson-TC
12cm



AG3-669
Adson-TC
12cm 1x2



AG3-670
Adson-TC
15cm 1x2



AG3-671
Adson
12cm



AG3-672
Adson
12cm 1x2



AG3-673
Adson
15cm 1x2



AG3-674
Adson Brown
12cm



AG3-675
Crane Kaplan
15cm



AG3-676
Crane Kaplan
15cm



AG3-677
Crane Kaplan
15cm



AG3-678
Hemingway
15cm



AG3-679
Miller
15cm Str.



AG3-680
Miller
15cm Cvd.



AG3-681
Dressing
Tweezer
14.5cm



AG3-682
College
Curved
Tip 15cm



AG3-683
Maryam
Tweezers
15cm



AG3-684
College
Tweezer
15cm



AG3-685
Perry
Tweezer
12.5cm



AG3-686
Self-locking
Tweezer
15cm

EXPLORERS


- Зонды
 - Sondes
 - Zahnsonden
- Зонди
 - Sondas
 - المستكشفين

Single Ended Probe Solid Handle 4.5mm Length 150mm



Explorers-Single-End


Double-Ended Explorers



EXPLORERS

- Зонды
 - Sondes
 - Zahnsonden
- Зонди
 - Sondas
 - المستكشفين

Probes Double Ended 6mm, length 170mm



Pointed Probes Double Ended Handle 6mm, length 170mm



EXPLORERS

- Зонды
 - Sondes
 - Zahnsonden
- Зонди
 - Sondas
 - المستكشفين

EXPLORERS INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM

PLUGGERS



AG3-724
5/7
21mm



AG3-725
9/11
21mm



AG3-726
1/2
Root Canal
Plugger



AG3-727
3/4
Root Canal
Plugger

.75

1.15

.60

1.10

ALLIED GROUP
of Industries



AG3-723
1/3
21mm



AG3-728
MA57
27mm



AG3-729
GP3
30mm



AG3-730
ISO
Size 30
24mm



AG3-731
ISO
Size 40
24mm



AG3-732
ISO
Size 50
24mm




AG3-733
ISO
Size 60
24mm

PERIODONTAL PROBES

- Зонды пародонтологические
 - Sondes paradontales
 - sonde parodontali
- зонди пародонтологічні
 - Sondas paradontales
 - المسابر اللثة

SPREADERS

HEAT CARRIERS



PERIODONTAL PROBES

- Зонды пародонтологические
 - Sondes parodontales
 - sonda parodontali
- зонди пародонтологічні
 - Sondas parodontales
 - المسابر اللثة

ENDODONTICS INSTRUMENTS

ANTERIOR



LUKS




SCALER



AG3-762
Gracey Curettes 1/2



AG3-763
Gracey Curettes 3/4




AG3-764
Gracey Curettes 5/6



AG3-765
Gracey Curettes 7/8



AG3-766
Gracey Curettes 9/10



AG3-767
Gracey Curettes 11/12



AG3-768
Gracey Curettes 13/14



AG3-769
Gracey Curettes 1/2



AG3-770
Gracey Curettes 3/4



AG3-771
Gracey Curettes 5/6



AG3-772
Gracey Curettes 7/8



AG3-773
Gracey Curettes 9/10



AG3-774
Gracey Curettes 11/12



AG3-775
Gracey Curettes 13/14



SCALER

GRACEY MINI FIVE

Hollow Handle 9.5mm, Length 170mm Ceramic Coated Handle



AG3-777
Fig. 3/4 Gracey 3M



AG3-778
Fig. 5/6 Gracey 3M



AG3-779
Fig. 7/8 Gracey 3M
Buccal and lingual/
Palatal surfaces



AG3-780
Fig. 11/12 Gracey 3M
Mesial



AG3-781
Fig. 13/14 Gracey 3M
Distal

AG3-776
Fig. 1/2
Gracey 3m



GRACEY AFTER FIVE

Hollow Handle 9.5mm, Length 170mm Ceramic Coated Handle



AG3-783
Gracey 3/4 AF



AG3-784
Gracey 5/6 AF



AG3-785
Gracey 7/8 AF



AG3-786
Gracey 11/12 AF



AG3-787
Gracey 13/14 AF

AG3-782
Gracey
1/2 AF



Gracey Curettes



$\frac{3...1}{3...1} | \frac{1...3}{1...3}$ **AG3-788** Fig. 1/2 GR - Gracey
Buccal surfaces



$\frac{5...1}{5...1} | \frac{1...5}{1...5}$ **AG3-789** Fig. 3/4 GR - Gracey
Lingual/Palatal surfaces



$\frac{5...1}{5...1} | \frac{1...5}{1...5}$ **AG3-790** Fig. 5/6 GR - Gracey



$\frac{8...4}{8...4} | \frac{4...8}{4...8}$ **AG3-791** Fig. 7/8 GR - Gracey
Buccal and lingual/
palatal surfaces



$\frac{8...4}{8...4} | \frac{4...8}{4...8}$ **AG3-792** Fig. 9/10 GR - Gracey
Buccal and lingual/
palatal surfaces more angled



$\frac{8...4}{8...4} | \frac{4...8}{4...8}$ **AG3-793** Fig. 11/12 GR - Gracey
Mesial




$\frac{8...4}{8...4} | \frac{4...8}{4...8}$ **AG3-794** Fig. 13/14 GR - Gracey
Distal



$\frac{8...4}{8...4} | \frac{4...8}{4...8}$ **AG3-795** Fig. 15/16 GR - Gracey
Mesial



$\frac{8...4}{8...4} | \frac{4...8}{4...8}$ **AG3-796** Fig. 17/18 GR - Gracey
Distal



SCALER

Sickle Scalars



AG3-797 Fig. U15/33 - Towner-Jaquette



AG3-798 Fig. U15/30 - Towner-Jaquette



AG3-799 Fig. H5/33



AG3-800 Fig. H6/H7



AG3-801 Fig. 30/33 - Jaquette



AG3-802 Fig. T2/T3 - Taylor



AG3-803 Fig. 204




AG3-804 Fig. 204S



AG3-805 Fig. 11/12 - MC - McCall
Mesial



AG3-806 Fig. 13/14 - MC - McCall
Distal





Straight and curved sickle scalars

Feature: 2 cutting edges curved facial component tipped ends



Universal curettes

are characterized by an angle of 90° between the terminal stem and the cutting component.



Gracey curettes

are characterized by an angle of 70° between the terminal stem and the cutting component.

Universal Curettes



AG3-807 Fig. 13S/14S - MC - McCall



AG3-808 Fig. 17S/18S - MC - McCall
Distal



AG3-809 Fig. M23



AG3-810 Fig. 204SD



AG3-811 Fig. YG7/8-Younger Good



AG3-812 Fig. L 1/2 - Langer



AG3-813 Fig. L 3/4 - Langer




AG3-814 Fig. L 5/6 - Langer



AG3-815 Fig. M 23A





AG3-816 Fig. M 23B




Straight and curved sickle scalers

Feature: 2 cutting edges curved facial component tipped ends



Universal curettes

are characterized by an angle of 90° between the terminal stem and the cutting component.



Gracey curettes

are characterized by an angle of 70° between the terminal stem and the cutting component.

Periodontal knives





0.9 MM

0.7 MM

Chisels
Diamond coated margin trimmers

Thanks to the diamond-coated hand instruments the margins are evenly levelled.

SCALER

Margin trimmers



AG3-835
1TH2 Hirt



AG3-836
5TH6 Hirt



AG3-837
3TH4 Hirt




AG3-838
7TH8 Hirt




AG3-839
62NS61

2.0 MM



AG3-840
37NS36

1.5 MM



AG3-841
66NS65

1.0 MM



AG3-842
79CI180 Black



AG3-843
77CI78 Black



AG3-844
279CI280 Black



AG3-845
277CI278 Black



AG3-846
179CI180 Black



AG3-847
177CI178 Black

CHISELS



AG3-848
1/2 Werdelstaedt



AG3-849
3/4 Werdelstaedt



AG3-850
7 Straight



AG3-851
5/6 Straight



AG3-852
7/10 Straight



AG3-853
8/9 Binangle



AG3-854
11/12 Binangle



AG3-855
40/41 Binangle

Tartar Remover Instruments Solid Handle Ø 6 mm



AG3-856
Fig. 1



AG3-857
Fig. 1L
Left



AG3-858
Fig. 1R
Right



AG3-859
Fig. 2L
left



AG3-860
Fig. 2R
Right



AG3-861



AG3-862



AG3-863

BURNISHERS CONSERVATIVE TREATMENT



FILLING INSTRUMENTS CONSERVATIVE TREATMENT



FILLING INSTRUMENTS CONSERVATIVE TREATMENT



BURNISHERS CONSERVATIVE TREATMENT

Ball Burnishers




Hatches




COMPOSITE INSTRUMENTS

Spatulas rounded and sharp




Spatulas Blade Thickness 0.22 mm



COMPOSITE INSTRUMENTS

Spatulas




Spatulas-Pluggers



EXCAVATOR INSTRUMENTS


Water drop Our excavators have particularly sharp blades, free of burrs and perfectly rounded.




FILLING INSTRUMENTS CONSERVATIVE TREATMENT



COMPOSITE BALL PLUGGERS



FILLING INSTRUMENTS CONSERVATIVE TREATMENT





FILLING INSTRUMENTS CONSERVATIVE TREATMENT



FILLING INSTRUMENTS CONSERVATIVE TREATMENT



FILLING INSTRUMENTS CONSERVATIVE TREATMENT



AG3-1005
WILLIAMS



AG3-1006
WILLIAMS



AG3-1007
WILLIAMS



AG3-1008



AG3-1009



AG3-1010



AG3-1011



AG3-1012




AG3-1013



AG3-1014

COMPOSITE PLUGGERS INSTRUMENTS





COMPOSITE PLUGGERS INSTRUMENTS



Serrated Pluggers



Ø 0.9 MM		Ø 1.2 MM
	AG3-1034 OP36M	
Ø 1.1 MM		Ø 1.6 MM
	AG3-1035 OP36D	
Ø 1.6 MM		Ø 2.2 MM
	AG3-1036 Op36	
Ø 1.0 MM		Ø 1.7 MM
	AG3-1037 Op52	
Ø 2.2 MM		Ø 2.6 MM
	AG3-1038 Zi51	
Ø 1.8 MM		Ø 2.2 MM
	AG3-1039 Op44	
Ø 1.7 MM		Ø 2.6 MM
	AG3-1040 1BR2	
Ø 2.6 MM		Ø 3.6 MM
	AG3-1041 2BR3	



COMPOSITE MODELLING INSTRUMENTS



AG3-1042 MF2C
Fuzzi



AG3-1043
Pj4



AG3-1044 MB3C
Besek



AG3-1045
HZ3C



AG3-1046
CPRO3C



AG3-1047 Mo1
Monari



AG3-1048 TH28C
Hirt



AG3-1049
HZ25C



AG3-1050
Ck10



AG3-1051
CK28BC



AG3-1052 Op37
Carver



AG3-1053
28BR29

EXCAVATORS

Round excavators

- Cariology**
- for the removal of softened, carious dentin
 - for the removal of temporary fillings
- Surgery**
- to sever the periodontal ligament prior to extraction
 - for the removal of granulation tissue after extraction
- Design: Round working end Curved cutting edge



Round Excavator, 1.0 mm
AG3-1054



Round Excavator, 1.5 mm
AG3-1055



Round Excavator, 2.0 mm
AG3-1056

Triple-angled back action excavators



Triple-angled Back Action Excavator, 1.0 mm
AG3-1057



Triple-angled Back Action Excavator, 1.5 mm
AG3-1058



Triple-angled Back Action Excavator, 2.0 mm
AG3-1059



Triple-angled Back Action Excavator, 3.5 mm
AG3-1060



FILLING INSTRUMENTS


- Инструменты для реставрации
 - Instruments pour plombage
 - Füllungsinstrumente (mit Kontrawinkel)
- Инструменти для реставрації
 - Instrumentos para obturaciones
 - ملء الصكوك

COMPOSITE/PLASTIC FILLING INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM

BLACK


PLUGGERS/CONDENSERS

Hollenback



MARQUETTE CONDENSER

MORTONSON



FILLING INSTRUMENTS

- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)
- Инструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

RESTORATIVE INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM Back Action Filling Instruments

Tanner Filling Instruments



AG3-1074
2T
Tanner

AG3-1075
3T
Tanner

AG3-1076
9/10
Back Action

AG3-1077
11/12
Back Action

AG3-1078
153
English Pattern
Serrated

AG3-1079
153
English Pattern
Serrated

BURNISHERS

BALL

BEAVERTAIL

ACORN

LADMORE



AG3-1073
1T
Tanner

AG3-1080
18

AG3-1081
26/27S

AG3-1082
27/29

AG3-1083
2

AG3-1084
2
Double-End

AG3-1085
21

AG3-1086
21B

AG3-1087
2
Ladmore

CARVERS INSTRUMENTS

- Инструменты для реставрации
- Instruments a modeler
- Modellierinstrumente
- Инструменти для реставрації
- Instrumentos para modelar
- ادوات تطعيم

P.K. THOMS WAXING INSTRUMENTS SOLID HANDLE 6mm




AG3-1095
Composite Instruments Blue
Titanium with Cassettes



AG3-1094
P.k. Thomas Waxing Instruments Kit

COMPOSITE INSTRUMENTS TITANIUM BLUE ENDODONTIC



CARVERS INSTRUMENTS

- Инструменты для реставрации
- Instruments a modeler
- Modellierinstrumente
- Инструменти для реставрації
- Instrumentos para modelar
- أدوات تقطیع



AG3-1102



AG3-1103



AG3-1104



AG3-1105



AG3-1106



Cassette For 5 Instruments
AG3-1107

FILLING INSTRUMENTS

- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)
- Инструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكونك

FILLING INSTRUMENTS (Gold Pointed) SOLID HANDLE 6MM, LENGTH 170MM



Filling instrument,
spatula 6 mm solid handle,
length 170 mm,

Filling instrument
6 mm
length 175 mm,


Filling instrument,
ball 6 mm
length 170 mm,

Filling instrument,
ball / spatula
6 mm
length 175 mm,

Filling instrument,
plan / plan 6 mm
length 165 mm,

COMPOSITE FILLING INSTRUMENTS

- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)
- Инструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك



AG3-1118

AG3-1119

AG3-1120

AG3-1121

AG3-1122

AG3-1123

AG3-1124

AG3-1125

AG3-1126

COMPOSITE FILLING INSTRUMENTS

- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)
- Инструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

Composite Gold Tip Instruments



MISCELLANEOUS PRODUCTS

- Прочие товары
 - Divers produits et services
 - Sonstige Produkte
- Інші товари
 - Varios Productos
 - المنتجات المتنوعة



AG3-1138
X Ray Film Hangers



AG3-1139
X Ray Film Hangers



AG3-1140
X Ray Film Hangers



AG3-1141
X Ray Film Hangers



AG3-1142
X Ray Film Hangers



AG3-1143
X Ray Film Hangers



AG3-1144
X Ray Film Hangers



AG3-1145



AG3-1146



AG3-1147



Dissecting Set
AG3-1148 Set of 9 Pcs.
AG3-1149 Set of 11 Pcs.



AG3-1150
Napkin Holder with Plastic Chain

AG3-1151
Napkin holder with metal chain



Mouth Props
AG3-1152
Small

Mouth Props
AG3-1153
Medium

Mouth Props
AG3-1154
Large

Key Rings



AG3-1155



AG3-1156



AG3-1157



AG3-1158



AG3-1159



AG3-1160

CEMENT SPATULAS

- Шпатели для цемента
 - Spatules à ciment
 - Zementspatel
- Шпателі для цементу
 - Espátulas para cemento
 - ملاعق الإسمنت



WAX KNIVES

- Ножи для воска
 - Couteaux a cire
 - Wachsmesser
- Ножи для воску
 - Cuchillos para cera
 - السكاكين الشمع

ALGINATE, PLASTER SPATULAS, WAX KNIVES PLASTER FOR KNIVES & PLASTER SHEAR



AG3-1183
Wooden Handle
Fig. 5
21.5 cm



AG3-1184
Wooden Handle
Fig. 7
19.5 cm curved



AG3-1185
Wooden Handle
Fig. 8
20 cm curved



AG3-1186
Wooden Handle
Fig. 9
22.5 cm curved



AG3-1187
Hook Knives



AG3-1188
Fig. 1
Lessmann
12.5 cm



AG3-1189
Fig. 2
Lessmann
17 cm



AG3-1190
Fig. 1
with small spoon
12.5 cm



AG3-1191
Fig. 2
with large spoon
17 cm



AG3-1192
Bacolite Handle
Wax Knives,
Small



AG3-1193
Bacolite Handle
Wax Knives,
Large



AG3-1194
Wooden Handle
Wax Knives,
Small



AG3-1195
Wooden Handle
Wax Knives,
Large



AG3-1196
Plaster Knife



AG3-1197
Fig. 1
17 cm



AG3-1198
Fig. 3
16 cm



AG3-1199
Fig. 4
16 cm



AG3-1200
Special Inox
19.5 cm

ARTICULATORS

- Артикуляторы
- Articular
- Artikulations

- Артикулятори
- Articuler
- آرتیکو لاتورس



AG3-1201

Technical Datas:

Intercondylar distance	120 mm.
condylar inclination	fixed at 30°
Bennet angle	fixed at 7°
Incisal table inclination	10°
Distance between mounting plates	85 mm.
Weight	400 gr.
Opening (item 5000)	125°
Opening (item 5010)	180°



AG3-1202

Technical Datas:

Intercondylar distance	110 mm.
condylar inclination	fixed at 30°
Bennet angle	fixed at 7°
Incisal table inclination	flat
Distance between mounting plates	85 mm.
Weight	450 gr.
Opening	180°



AG3-1203

Stainless Steel Finishing

AG3-1204

Royal Blue



AG3-1205



AG3-1206



AG3-1207



AG3-1208



AG3-1209



AG3-1210



AG3-1211



AG3-1212




AG3-1213



AG3-1214
Nakawalni
Dental Laboratory
Instrument

DENTAL PHOTOGRAPHIC MIRRORS



AG3-1215 Adult Palatal /
AG3-1216 Child Palatal

AG3-1217
Buccal

AG3-1218 Adult Palatal /
AG3-1219 Child Palatal

AG3-1220 Adult Palatal /
AG3-1221 Adult Buccal




AG3-1222 Lingual upper /
AG3-1223 Lingual Lower

AG3-1224 Adult Buccal /
AG3-1225 Child Buccal

AG3-1226 Child Palatal /
AG3-1227 Child Buccal

AG3-1228
Palatal-Adult



AG3-1229
Palatal-Adult Large

AG3-1230
Palatal-Child

AG3-1231
Buccal-Small

AG3-1232
Buccal

AG3-1233
Lingual

AG3-1234
Lingual-Large